



1. THE SLIDING ARM CAN BE POSITIONED ON THE CONNECTING ROD AS REQUIRED WITHOUT CHANGING THE DIRECTION OF THE COAXIAL PLATE;
2. LARGER OPTICAL COMPONENTS CAN BE MOUNTED OUTSIDE THE COAXIAL SYSTEM USING LONGER MOUNTING ADAPT

DRAWING PROJECTION			CSM2-CL		Zolix			
	NAME	DATE	60mm CYLINDRICAL LENS MOUNT		MATERIAL	WEIGHT	SCALE	REV
DRAWN	NIC	4th/August/2025			FOLLOW THE PARTS	0.093	1:1	V1
APPROVAL	DICK	4th/August/2025	THE DRAWINGS ARE ONLY FOR TYPE SELECTION INFORMATION		shop.zolix.com.cn			